



# New SuperLine II

## Product Catalog

**Dentium**



# CONTENTS

## Introduction

S.L.A. Surface	04
SuperLine Characteristics	06
SuperLine Color Coding by Diameter	08
SuperLine Fixture	09

## Surgical Components

Cover Screw	11
Healing Abutment	12
Scan Body	13

## Prosthetic Procedure 1

Dual / Combi Abutment - Abutment Level Impression	14
Combi Abutment	15
Dual Abutment[Hex]	16
Dual Abutment[Non-hex]	17
Abutment Level Impression Components	18
Restorative Kit	19

## Prosthetic Procedure 2

Dual / Custom / Dual Milling / Angled / Direct-Casting / Metal-Casting / Temporary (Plastic & Ti) Abutment - Fixture Level Impression	20
Fixture Level Impression Components	21
Custom Abutment	23
Titanium Abutment Blank	24
Dual Milling Abutment[Ti-G4]	25
Dual Milling Abutment[Ti-G4]	26
Angled Abutment[Ti-G4 / 15°]	27
Angled Abutment[Ti-G4 / 25°]	28
Metal-Casting Abutment	29
Temporary Abutment	30

## Prosthetic Procedure 3

Screw Abutment - Abutment Level Impression	31
Screw Abutment	32
Angled Screw Abutment	33
Screw Abutment Impression Components	34


## Prosthetic Procedure 4

Overdenture Procedure	
- Positioner / Mini Ball / Magnetic Attachment	38
Positioner	39
Mini Ball Attachment	41
Magnetic Attachment (Dome Type)	42
Magnetic Attachment (Flat Type)	43

## Instruments

Surgical Kit (UXIF)	44
Surgical Kit (USNF)	45
Surgical Kit (UXSF)	46
Digital Guide Full Kit	47
Surgical Kit (UXIFN)	48
Surgical Kit (XSIF)	49
Drill Stopper Kit	50
Drill	51
Stopper	55
Instruments	56
Dentium Advanced Sinus Kit (DASK)	59
Sinus Bur Kit / Sinus Kit	60
Dask / Sinus Bur Kit / Sinus Kit	61
Dask Simple	62
Ridge Spreader Kit	65
Osteotome Kit	66
Trephine Kit	67
Prosthetic Kit	68
Planning Kit	69
Prosthetic and Laboratory Instrument	70
Help Kit	72

<b>Surgical Manual</b>	73
<b>Prosthesis Manual</b>	83







## S.L.A. Surface

**S.L.A.** (Sandblasting with large grit and acid etching)

- Higher bone-to-implant contact.
- Faster bone formation on the surface.

*In vivo test*







***Super*** Líne



# SuperLine Characteristics

## "Immediate Implantation with Excellent Bone Response"

- Higher stabilization in extraction socket
- Early loading in upper posterior
- Harmony with anatomy
- Sharp & fast insertion



### Selection Guideline


Ideal emergence profile for each tooth





# SuperLine Characteristics

## A Double-threaded Design & Extended Cutting Edge




- Improved self-tapping ability
- Allows more control over the depth of fixture placement
- Alleviates the occurrence of over-torque during placement in dense bone
- Sharpened thread design promotes better initial stability in soft bone
- Easy & fast insertion can be done due to double threaded straight body design

## B Biological Connection




- The conical hex connection between implant and abutment interface ensures hermetic sealing.
- The biologic connection distributes the load to the fixture evenly. Therefore it helps minimize micro-movement and marginal bone loss.
- All implant diameters share the same internal hex.

## C Tapered Design




- Tapered design may harmonize with surrounding bone anatomically.
- An ideal design for greater stability during sinus surgeries.

## D Osseointegration



- The greater distance between the threads may promote early osseointegration

## E Prosthesis










- One abutment screw fits all abutments and fixture platforms.
- Single abutment connection is used for all implant diameters.
- One hex screw driver fits all abutment screws



# SuperLine Color Coding by Diameter

## Color Coding by Diameter

• Cover screw is not included. (Unit: mm)

Platform		Cap Color						
			Yellow	Green	Blue	Red	Orange	Violet
Fixture SuperLine (Mount Free)								
	PlatformØ Fixture Platform Diameter		3.6	4.0	4.5	5.0	6.0	7.0
	BodyØ Fixture Body Diameter		3.4	3.8	4.3	4.8	5.0	5.8
	L : 7 Fixture Bevel Height		1.5	1.5	1.5	1.5	1.5	1.5
	L : 8, 10, 12, 14 Fixture Bevel Height		0	0.1	0.3	0.4	0.7	1.0





# SuperLine Fixture

Unit: mm, Scale 1 : 1.5 / mm



**Platform Ø3.6 | Body Ø3.4**

L	Art. No.
7	FXS 36 <b>07</b>
8	FXS 36 <b>08</b>
10	FXS 36 <b>10</b>
12	FXS 36 <b>12</b>
14	FXS 36 <b>14</b>




**Platform Ø4.0 | Body Ø3.8**

L	Art. No.
7	FXS 40 <b>07</b>
8	FXS 40 <b>08</b>
10	FXS 40 <b>10</b>
12	FXS 40 <b>12</b>
14	FXS 40 <b>14</b>




**Platform Ø4.5 | Body Ø4.3**

L	Art. No.
7	FXS 45 <b>07</b>
8	FXS 45 <b>08</b>
10	FXS 45 <b>10</b>
12	FXS 45 <b>12</b>
14	FXS 45 <b>14</b>



**Platform Ø5.0 | Body Ø4.8**

L	Art. No.
7	FXS 50 <b>07</b>
8	FXS 50 <b>08</b>
10	FXS 50 <b>10</b>
12	FXS 50 <b>12</b>
14	FXS 50 <b>14</b>




※Note: 1) To prevent damage to the Implant driver or fixture, do not over torque during fixture insertion

# SuperLine Fixture

Unit: mm, Scale 1 : 1.5 / mm


Platform Ø6.0 | Body Ø5.0

L	Art. No.
7	FXS 60 <b>07</b>
8	FXS 60 <b>08</b>
10	FXS 60 <b>10</b>
12	FXS 60 <b>12</b>



Platform Ø7.0 | Body Ø5.8

L	Art. No.
7	FXS 70 <b>07</b>
8	FXS 70 <b>08</b>
10	FXS 70 <b>10</b>
12	FXS 70 <b>12</b>






# Cover Screw

- Single use only

Unit: mm, Scale 1 : 1.5 / mm



CS36 and FXS4510

Color: Green


Art. No.	CS36
----------	------



# Healing Abutment

- Single use only

Unit: mm, Scale 1 : 1.5 / mm




HAB453050L and FXS4510

※ Hex driver: Use no more than 5N·cm torque when screwing a cover screw to a fixture.  
If hex is worn, slot on the head of the product can be used to rotate it.

# Healing Abutment


## Diameter Ø4.0

G/H	H	Art. No.
2.0	2.0	HAB 40 <b>20</b> 20 L
2.0	3.5	HAB 40 <b>20</b> 35 L
3.0	5.0	HAB 40 <b>30</b> 50 L
4.0	7.0	HAB 40 <b>40</b> 70 L




## Diameter Ø4.5

G/H	H	Art. No.
2.0	2.0	HAB 45 <b>20</b> 20 L
2.0	3.5	HAB 45 <b>20</b> 35 L
3.0	5.0	HAB 45 <b>30</b> 50 L
4.0	7.0	HAB 45 <b>40</b> 70 L




## Diameter Ø5.5

G/H	H	Art. No.
2.0	2.0	HAB 55 <b>20</b> 20 L
2.0	3.5	HAB 55 <b>20</b> 35 L
3.0	5.0	HAB 55 <b>30</b> 50 L
4.0	7.0	HAB 55 <b>40</b> 70 L



## Diameter Ø6.5

G/H	H	Art. No.
2.0	2.0	HAB 65 <b>20</b> 20 L
2.0	3.5	HAB 65 <b>20</b> 35 L
3.0	5.0	HAB 65 <b>30</b> 50 L
4.0	7.0	HAB 65 <b>40</b> 70 L



## Diameter Ø7.5 / 8.5 / 9.5

G/H	H	Art. No.
3.0	5.0	HAB 75 <b>30</b> 50 L
3.0	5.0	HAB 85 <b>30</b> 50 L
3.0	5.0	HAB 95 <b>30</b> 50 L



※ Hex driver: Use no more than 10N·cm of torque when screwing a healing abutment to a fixture.  
If hex is worn, slot on the head of the product can be used to rotate it.



# Scan Body



- Single use only, provided sterile
- Tissue contour design
- Abutment screw included
- Material: Ti-6Al-4V ELI
- **NOTE:** Use no more than 10 N·cm of torque when tightening the Scan body

Unit: mm    Scale 1.5:1


Diameter	G/H	Height	Type	Ref. No.
Ø 4.0	2.5	6.0	Hex	IHAB 40 06 TH
	2.5	6.0	Non-hex	IHAB 40 06 TN
	3.0	8.0	Hex	IHAB 40 08 TH
	3.0	8.0	Non-hex	IHAB 40 08 TN
	3.5	10.0	Hex	IHAB 40 10 TH
	3.5	10.0	Non-hex	IHAB 40 10 TN
Ø 5.0	2.5	6.0	Hex	IHAB 50 06 TH
	2.5	6.0	Non-hex	IHAB 50 06 TN
	3.0	8.0	Hex	IHAB 50 08 TH
	3.0	8.0	Non-hex	IHAB 50 08 TN
	3.5	10.0	Hex	IHAB 50 10 TH
	3.5	10.0	Non-hex	IHAB 50 10 TN
Ø 6.0	2.5	6.0	Hex	IHAB 60 06 TH
	2.5	6.0	Non-hex	IHAB 60 06 TN
	3.0	8.0	Hex	IHAB 60 08 TH
	3.0	8.0	Non-hex	IHAB 60 08 TN
	3.5	10.0	Hex	IHAB 60 10 TH
	3.5	10.0	Non-hex	IHAB 60 10 TN



# Prosthetic Procedure 1

Impression Technique and Restoration Selection


## Dual / Combi Abutment






# Combi Abutment

Unit: mm, Scale 1 : 1 / mm




CAB5535SL and FXS4510



CAB5535L and FXS4510


## Diameter Ø4.5

G/H	Type	Art. No.	Type	Art. No
1.0	Short	CAB 45 <b>10</b> SL	Long	CAB 45 <b>10</b> L
1.5		CAB 45 <b>15</b> SL		CAB 45 <b>15</b> L
2.5		CAB 45 <b>25</b> SL		CAB 45 <b>25</b> L
3.5		CAB 45 <b>35</b> SL		CAB 45 <b>35</b> L
4.5		CAB 45 <b>45</b> SL		CAB 45 <b>45</b> L
5.5		CAB 45 <b>55</b> SL		CAB 45 <b>55</b> L




## Diameter Ø5.5

G/H	Type	Art. No.	Type	Art. No
1.5	Short	CAB 55 <b>15</b> SL	Long	CAB 55 <b>15</b> L
2.5		CAB 55 <b>25</b> SL		CAB 55 <b>25</b> L
3.5		CAB 55 <b>35</b> SL		CAB 55 <b>35</b> L
4.5		CAB 55 <b>45</b> SL		CAB 55 <b>45</b> L
5.5		CAB 55 <b>55</b> SL		CAB 55 <b>55</b> L



## Diameter Ø6.5

G/H	Type	Art. No.	Type	Art. No
1.5	Short	CAB 65 <b>15</b> SL	Long	CAB 65 <b>15</b> L
2.5		CAB 65 <b>25</b> SL		CAB 65 <b>25</b> L
3.5		CAB 65 <b>35</b> SL		CAB 65 <b>35</b> L
4.5		CAB 65 <b>45</b> SL		CAB 65 <b>45</b> L
5.5		CAB 65 <b>55</b> SL		CAB 65 <b>55</b> L




※Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the combi abutment with fixture.

# Dual Abutment [Hex]

- Abutment screw is included.

Unit: mm, Scale 1 : 1.5 / mm



DAB5535HL and FXS4510

## Diameter Ø4.5 | Hex

G/H	Art. No.
1.0	DAB 45 <b>10</b> HL
1.5	DAB 45 <b>15</b> HL
2.5	DAB 45 <b>25</b> HL
3.5	DAB 45 <b>35</b> HL
4.5	DAB 45 <b>45</b> HL
5.5	DAB 45 <b>55</b> HL



## Diameter Ø5.5 | Hex

G/H	Art. No.
1.5	DAB 55 <b>15</b> HL
2.5	DAB 55 <b>25</b> HL
3.5	DAB 55 <b>35</b> HL
4.5	DAB 55 <b>45</b> HL
5.5	DAB 55 <b>55</b> HL



## Diameter Ø6.5 | Hex

G/H	Art. No.
1.5	DAB 65 <b>15</b> HL
2.5	DAB 65 <b>25</b> HL
3.5	DAB 65 <b>35</b> HL
4.5	DAB 65 <b>45</b> HL
5.5	DAB 65 <b>55</b> HL




※Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the Dual abutment with fixture.



# Dual Abutment [Non-hex ]

- Abutment screw is included.


Unit: mm, Scale 1 : 1.5 / mm



DAB5535NL and FXS4510


## Diameter Ø4.5 | Non-hex

G/H	Art. No.
1.0	DAB 45 <b>10</b> NL
1.5	DAB 45 <b>15</b> NL
2.5	DAB 45 <b>25</b> NL
3.5	DAB 45 <b>35</b> NL
4.5	DAB 45 <b>45</b> NL
5.5	DAB 45 <b>55</b> NL




## Diameter Ø5.5 | Non-hex

G/H	Art. No.
1.5	DAB 55 <b>15</b> NL
2.5	DAB 55 <b>25</b> NL
3.5	DAB 55 <b>35</b> NL
4.5	DAB 55 <b>45</b> NL
5.5	DAB 55 <b>55</b> NL



## Diameter Ø6.5 | Non-hex

G/H	Art. No.
1.5	DAB 65 <b>15</b> NL
2.5	DAB 65 <b>25</b> NL
3.5	DAB 65 <b>35</b> NL
4.5	DAB 65 <b>45</b> NL
5.5	DAB 65 <b>55</b> NL



※Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the Dual abutment with fixture.

# Abutment Level Impression Components

Unit: mm, Scale 1 : 1 / mm

## Comfort Cap | Snap on

Type	Diameter	Art. No.
Short	Ø4.5	CCC 45 CS
	Ø5.5	CCC 55 CS
	Ø6.5	CCC 65 CS
Long	Ø4.5	CCC 45 C
	Ø5.5	CCC 55 C
	Ø6.5	CCC 65 C




## Impression Coping

Diameter	Art. No.
Ø4.5	CIC 45 L
Ø5.5	CIC 55 L
Ø6.5	CIC 65 L




## Lab Analog

Type	Diameter	Art. No.
Short	Ø4.5	CAN 45 SL
	Ø5.5	CAN 55 SL
	Ø6.5	CAN 65 SL
Long	Ø4.5	CAN 45 LL
	Ø5.5	CAN 55 LL
	Ø6.5	CAN 65 LL



## Burn-out Cylinder

Type	Diameter	Art. No.
Single	Ø4.5	CBC 45 SL
	Ø5.5	CBC 55 SL
	Ø6.5	CBC 65 SL
Bridge	Ø4.5	CBC 45 BL
	Ø5.5	CBC 55 BL
	Ø6.5	CBC 65 BL



# Restorative Kit



## Combi & Dual Abutment


Art. No	Lab. Components				
	Comfort Cap	Impression Coping	Analog	Burn-out Cylinder	
XSDAB <b>45 S</b> XSDAB <b>45</b>	CCC 45 CS CCC 45 C	CIC 45 L	CAN 45 SL CAN 45 LL	CBC 45 SL	CBC 45 BL
XSDAB <b>55 S</b> XSDAB <b>55</b>	CCC 55 CS CCC 55 C	CIC 55 L	CAN 55 SL CAN 55 LL	CBC 55 SL	CBC 55 BL
XSDAB <b>65 S</b> XSDAB <b>65</b>	CCC 65 CS CCC 65 C	CIC 65 L	CAN 65 SL CAN 65 LL	CBC 65 SL	CBC 65 BL



# Prosthetic Procedure 2

Impression Technique and Restoration Selection

Dual / Custom / Dual Milling / Angled /  
Metal-Casting / Temporary (Plastic & Ti) Abutment




# Fixture Level Impression Components

• Impression coping screw is included with Impression coping.

Unit: mm, Scale 1 : 1.5 / mm


## Impression Coping Pick-up Ø4.0

Size	Type	Art. No.
Short	<b>Hex</b>	DPU 40 11 HL
Short	<b>Non-hex</b>	DPU 40 11 NL
Long	<b>Hex</b>	DPU 40 15 HL
Long	<b>Non-hex</b>	DPU 40 15 NL




## Impression Coping Pick-up Ø4.5

Size	Type	Art. No.
Short	<b>Hex</b>	DPU 45 11 HL
Short	<b>Non-hex</b>	DPU 45 11 NL
Long	<b>Hex</b>	DPU 45 15 HL
Long	<b>Non-hex</b>	DPU 45 15 NL




## Impression Coping Pick-up Ø5.5

Size	Type	Art. No.
Short	<b>Hex</b>	DPU 55 11 HL
Short	<b>Non-hex</b>	DPU 55 11 NL
Long	<b>Hex</b>	DPU 55 15 HL
Long	<b>Non-hex</b>	DPU 55 15 NL




## Impression Coping Pick-up Ø6.5

Size	Type	Art. No.
Short	<b>Hex</b>	DPU 65 11 HL
Short	<b>Non-hex</b>	DPU 65 11 NL
Long	<b>Hex</b>	DPU 65 15 HL
Long	<b>Non-hex</b>	DPU 65 15 NL



## Impression Coping Pick-up Screw

Size	Art. No.
Short	DPS 11
Long	DPS 15



## Analogue

Application (BodyØ)	Art. No.
<b>3.4 / Ø3.8 / Ø4.3 / Ø4.8 / Ø5.8</b>	DAN38




# Fixture Level Impression Components

• Impression coping screw is included with Impression coping. Unit: mm, Scale 1 : 1.5 / mm


## Impression Coping Transfer Ø4.0

Size	Type	Art. No.
Short	<b>Hex</b>	DTF <b>40</b> 11 <b>HL</b>
Short	<b>Non-hex</b>	DTF <b>40</b> 11 <b>NL</b>
Long	<b>Hex</b>	DTF <b>40</b> 15 <b>HL</b>
Long	<b>Non-hex</b>	DTF <b>40</b> 15 <b>NL</b>




## Impression Coping Transfer Ø4.5

Size	Type	Art. No.
Short	<b>Hex</b>	DTF <b>45</b> 11 <b>HL</b>
Short	<b>Non-hex</b>	DTF <b>45</b> 11 <b>NL</b>
Long	<b>Hex</b>	DTF <b>45</b> 15 <b>HL</b>
Long	<b>Non-hex</b>	DTF <b>45</b> 15 <b>NL</b>




## Impression Coping Transfer Ø5.5

Size	Type	Art. No.
Short	<b>Hex</b>	DTF <b>55</b> 11 <b>HL</b>
Short	<b>Non-hex</b>	DTF <b>55</b> 11 <b>NL</b>
Long	<b>Hex</b>	DTF <b>55</b> 15 <b>HL</b>
Long	<b>Non-hex</b>	DTF <b>55</b> 15 <b>NL</b>




## Impression Coping Transfer Ø6.5

Size	Type	Art. No.
Short	<b>Hex</b>	DTF <b>65</b> 11 <b>HL</b>
Short	<b>Non-hex</b>	DTF <b>65</b> 11 <b>NL</b>
Long	<b>Hex</b>	DTF <b>65</b> 15 <b>HL</b>
Long	<b>Non-hex</b>	DTF <b>65</b> 15 <b>NL</b>



## Impression Coping Transfer Screw

Size	Art. No.
Short	DTS 11
Long	DTS 15



## Analog I for SuperLine and Implantium

Application (BodyØ)	Art. No.
<b>3.4 / Ø3.8 / Ø4.3 / Ø4.8 / Ø5.8</b>	DAN38





# Custom Abutment


- Abutment screw is included.

Unit: mm, Scale 1 : 1.5 / mm




## Diameter $\varnothing 4.5$

G/H	Type	Art. No.
0.5	Hex	CDAB 45 05 <b>H</b>
0.5	Non-hex	CDAB 45 05 <b>N</b>
1.5	Hex	CDAB 45 15 <b>H</b>
1.5	Non-hex	CDAB 45 15 <b>N</b>



## Diameter $\varnothing 5.5$

G/H	Type	Art. No.
1.0	Hex	CDAB 55 10 <b>H</b>
1.0	Non-hex	CDAB 55 10 <b>N</b>
2.0	Hex	CDAB 55 20 <b>H</b>
2.0	Non-hex	CDAB 55 20 <b>N</b>



※Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the custom abutment with fixture.

# Titanium Abutment Blanks

Unit: mm, Scale 1 : 0.8

1




- Two abutment screws included

Material : Ti-6Al-4V ELI


Diameter	Type	Ref. No.
Ø 10.0	Hex	CMAB 10 HAL
	Non-hex	CMAB 10 NA



Diameter	Type	Ref. No.
Ø 14.0	Hex	CMAB 14 HAL
	Non-hex	CMAB 14 NA




Diameter	Type	Ref. No.
Ø 15.0	Hex	CMAB 15 HAL
	Non-hex	CMAB 15 NA



# Dual Milling Abutment [Ti-G4 ]


• Abutment screw is included.

Unit:mm, Scale 1: 1.5 / mm




## Diameter $\varnothing 4.0$

G/H	Type	Art. No.
1.0	Hex	DAB 40 105 HL
1.0	Non-hex	DAB 40 105 NL




## Diameter $\varnothing 4.5$

G/H	Type	Art. No.
1.5	Hex	DAB 45 156 HL
1.5	Non-hex	DAB 45 156 NL



## Diameter $\varnothing 5.5$

G/H	Type	Art. No.
1.5	Hex	DAB 55 157 HL
1.5	Non-hex	DAB 55 157 NL
2.5	Hex	DAB 55 257 HL
2.5	Non-hex	DAB 55 257 NL



※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the milling abutment with fixture.




# Dual Milling Abutment [Ti-G4 ]

• Abutment screw is included.

Unit:mm, Scale 1: 1.5 / mm


## Diameter Ø6.5

G/H	Type	Art. No.
1.5	Hex	DAB 65 158 <b>HL</b>
1.5	Non-hex	DAB 65 158 <b>NL</b>
2.5	Hex	DAB 65 258 <b>HL</b>
2.5	Non-hex	DAB 65 258 <b>NL</b>
3.5	Hex	DAB 65 358 <b>HL</b>
3.5	Non-hex	DAB 65 358 <b>NL</b>



## Diameter Ø7.5

G/H	Type	Art. No.
2.5	Hex	DAB 75 259 <b>HL</b>
2.5	Non-hex	DAB 75 259 <b>NL</b>
3.5	Hex	DAB 75 359 <b>HL</b>
3.5	Non-hex	DAB 75 359 <b>NL</b>




※Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the milling abutment with fixture.

# Angled Abutment [Ti-G4 / 15° ]

• Abutment screw is included.


Unit: mm, Scale 1 : 1.5 / mm



AAB154515HL and FXS4510


## Diameter Ø4.5 | Angled 15°

G/H	Type	Art. No.
1.5	Hex	AAB 15 45 15 HL
1.5	Non-hex	AAB 15 45 15 NL
2.5	Hex	AAB 15 45 25 HL
2.5	Non-hex	AAB 15 45 25 NL
3.5	Hex	AAB 15 45 35 HL
3.5	Non-hex	AAB 15 45 35 NL



## Diameter Ø5.5 | Angled 15°

G/H	Type	Art. No.
1.5	Hex	AAB 15 55 15 HL
1.5	Non-hex	AAB 15 55 15 NL
2.5	Hex	AAB 15 55 25 HL
2.5	Non-hex	AAB 15 55 25 NL
3.5	Hex	AAB 15 55 35 HL
3.5	Non-hex	AAB 15 55 35 NL




※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the angled abutment with fixture.

# Angled Abutment [Ti-G4 / 25° ]

• Abutment screw is included.


Unit: mm, Scale 1 : 1.5 / mm



AAB254515HL and FXS4510


Diameter Ø4.5 | Angled 25°

G/H	Type	Art. No.
1.5	Hex	AAB 25 45 15 HL
1.5	Non-hex	AAB 25 45 15 NL
2.5	Hex	AAB 25 45 25 HL
2.5	Non-hex	AAB 25 45 25 NL
3.5	Hex	AAB 25 45 35 HL
3.5	Non-hex	AAB 25 45 35 NL



Diameter Ø5.5 | Angled 25°

G/H	Type	Art. No.
1.5	Hex	AAB 25 55 15 HL
1.5	Non-hex	AAB 25 55 15 NL
2.5	Hex	AAB 25 55 25 HL
2.5	Non-hex	AAB 25 55 25 NL
3.5	Hex	AAB 25 55 35 HL
3.5	Non-hex	AAB 25 55 35 NL



※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the angled abutment with fixture.




# Metal-Casting Abutment



RAB45CH and FXS4510

Metal-Casting Abutment | Co-Cr

G/H	Type	Art. No.
1.0	Hex	RAB 45 CH
1.0	Non-hex	RAB 45 CN




※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the meta casting abutment with fixture.

# Temporary Abutment


• Abutment screw is included.

Unit: mm, Scale 1 : 1.5 / mm




## Ti-Temporary Abutment

Diameter	G/H	Type	Art. No.
Ø4.5	1.0	Hex	RAB 45 THG
Ø4.5	1.0	Non-hex	RAB 45 TNG



## Plastic Temporary Abutment

Diameter	G/H	Type	Art. No.
Ø4.5	2.0	hex	RAB 45 20 PHL
Ø4.5	2.0	Non-hex	RAB 45 20 PNL
Ø5.5	2.0	hex	RAB 55 20 PHL
Ø5.5	2.0	Non-hex	RAB 55 20 PNL
Ø6.5	2.0	hex	RAB 65 20 PHL
Ø6.5	2.0	Non-hex	RAB 65 20 PNL




※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the temporary abutment with fixture.

# Prosthetic Procedure 3

Impression Technique and Restoration Selection

## Screw Abutment




# Screw Abutment

Unit: mm, Scale 1 : 1.5 / mm




Delivery Holder



SAB5535L and FXS4510


## Diameter Ø4.5

G/H	Art. No.
1.0	SAB45 10 L
1.5	SAB45 15 L
2.5	SAB45 25 L
3.5	SAB45 35 L
4.5	SAB45 45 L
5.5	SAB45 55 L



## Diameter Ø5.5

G/H	Art. No.
1.5	SAB55 15 L
2.5	SAB55 25 L
3.5	SAB55 35 L
4.5	SAB55 45 L
5.5	SAB55 55 L




※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the screw abutment with fixture.




# Angled Screw Abutment

Unit: mm, Scale 1 : 1.5 / mm




Ti-Delivery Holder

ASAH
------




Angled Screw Abutment Ø4.5

G/H	Type	Art. No.
1.0	Hex	ASA 45 15 10 18 H
1.0	Non-hex	ASA 45 15 10 18 N
1.0	Hex	ASA 45 30 10 18 H
1.0	Non-hex	ASA 45 30 10 18 N
2.0	Hex	ASA 45 15 20 18 H
2.0	Non-hex	ASA 45 15 20 18 N
2.0	Hex	ASA 45 30 20 18 H
2.0	Non-hex	ASA 45 30 20 18 N




Angled Screw Abutment Ø5.5

G/H	Type	Art. No.
1.5	Hex	ASA 55 15 15 18 H
1.5	Non-hex	ASA 55 15 15 18 N
1.5	Hex	ASA 55 30 15 18 H
1.5	Non-hex	ASA 55 30 15 18 N
3.0	Hex	ASA 55 15 30 18 H
3.0	Non-hex	ASA 55 15 30 18 N
3.0	Hex	ASA 55 30 30 18 H
3.0	Non-hex	ASA 55 30 30 18 N



Angled Screw Abutment Screw

ASASC 20 23
-------------





※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the angled screw abutment with fixture.

# Screw Abutment Impression Components

Unit: mm, Scale 1 : 1.5 / mm



## Impression Coping Pick-up

Diameter	Type		Art. No.
Ø4.5	Single	Hex	SPU 45 SL
Ø4.5	Bridge	Non-hex	SPU 45 BL
Ø5.5	Single	Hex	SPU 55 SL
Ø5.5	Bridge	Non-hex	SPU 55 BL




## Impression Coping Transfer

Diameter	Type		Art. No.
Ø4.5	Single	Hex	STF 45 SL
Ø4.5	Bridge	Non-hex	STF 45 BL
Ø5.5	Single	Hex	STF 55 SL
Ø5.5	Bridge	Non-hex	STF 55 BL



## Impression Coping Screw

Type	Art. No.
Pick-up	SPS 09
Transfer	STS 09



# Screw Abutment Components

## Scan Comfort Cap | for Screw Abutment

Diameter	Length	Type	Ref. No.
Ø 4.5	5.0	Hex Non-hex	ISC 45 ST ISC 45 BT
Ø 5.5		Hex Non-hex	ISC 55 ST ISC 55 BT

• Ti-retaining screw included



## Scan Comfort Cap | for Screw Abutment


Diameter	Length	Type	Ref. No.
Ø 4.5	5.0	Hex Non-hex	ISC 45 SP ISC 45 BP
Ø 5.5		Hex Non-hex	ISC 55 SP ISC 55 BP

• Ti-retaining screw (SRS 18T) included



## Analogue

Diameter	Art. No.
Ø4.5	SAN 45 L
Ø5.5	SAN 55 L



## Digital Analogue | for Screw Abutment

Diameter	Art. No.
Ø 4.5	SAN 45 D
Ø 5.5	SAN 55 D



## Polishing Protector

Diameter	Art. No.
Ø4.5	SPP 45 L
Ø5.5	SPP 55 L



## Ti-Retaining Screw


SRS 18 T
----------



# Screw Abutment Impression Components

## Temporary Cylinder | Ti-Cylinder


Diameter	Type		Art. No.
Ø4.5	Single	Hex	STC 45 SG
Ø4.5	Bridge	Non-hex	STC 45 BG
Ø5.5	Single	Hex	STC 55 SG
Ø5.5	Bridge	Non-hex	STC 55 BG



### Ti-CylinderMaterial : Ti-6AL-4V ELI

Diameter	Type	Art. No.
Ø 4.5	Hex	STA 45 S
	Non-hex	STA 45 B
Ø 5.5	Hex	STA 55 S
	Non-hex	STA 55 B

• Ti-retaining screw included




# Screw Abutment Impression Components

Unit: mm, Scale 1 : 1.5 / mm


## Metal Cylinder | Co-Cr

Diameter	Type		Art. No.
Ø4.5	Single	Hex	SGC 45 CSL
Ø4.5	Bridge	Non-hex	SGC 45 CBL
Ø5.5	Single	Hex	SGC 55 CSL
Ø5.5	Bridge	Non-hex	SGC 55 CBL




## Burn-out Cylinder

Diameter	Type		Art. No.
Ø4.5	Single	Hex	SBC 45 SL
Ø4.5	Bridge	Non-hex	SBC 45 BL
Ø5.5	Single	Hex	SBC 55 SL
Ø5.5	Bridge	Non-hex	SBC 55 BL



## Digital Analog

Implant platform Ø	Ref. No.
All diameters	DAN 38 D





# Prosthetic Procedure 4

Impression Technique and Restoration Type


## Overdenture Procedure

### Positioner / Mini Ball / Magnetic Attachment



# Positioner Abutment


Unit: mm, Scale 1 : 1.5 / mm



FSMH and PAB3520 and FXS4510

## Positioner Abutment

G/H	Art. No.
1.0	PAB 35 <b>10</b>
2.0	PAB 35 <b>20</b>
3.0	PAB 35 <b>30</b>
4.0	PAB 35 <b>40</b>
5.0	PAB 35 <b>50</b>
6.0	PAB 35 <b>60</b>
7.0	PAB 35 <b>70</b>



## Positioner Impression Coping

PIC
-----



## Positioner Analog

PAN
-----



※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the positioner abutment with fixture.

# Positioner

Unit: mm, Scale 1 : 1.5 / mm

## Positioner Socket Set

Art. No.	FSMHS(Tilting Type $\pm 10^{\circ}$ )
	FSMHSN(Non Tilting Type $\pm 5^{\circ}$ )




## Positioner Metal Socket

Art. No.	FSMH
----------	------



## Positioner Plastic Socket

Application	Art. No.
Tilting Type $\pm 10^{\circ}$	MSHP (Blue)
	MSMP (Orange)
	MSLP (Ivory)
	MSOP (White)
Non Tilting Type $\pm 5^{\circ}$	MSHPN (Blue)
	MSMPN (Orange)
	MSLPN (Ivory)
	MSOP (White)



## Positioner Block Out Spacer

Art. No.	PBOS
----------	------



## Positioner Core Tool

Art. No.	XPCT
----------	------



# Mini Ball Attachment


Unit: mm, Scale 1 : 1.5 / mm



BPF3 and BAB352018 and FXS4510


## Mini Ball Abutment

G/H	Art. No.
0	BAB 35 <b>00</b> 18
1.0	BAB 35 <b>10</b> 18
2.0	BAB 35 <b>20</b> 18
3.0	BAB 35 <b>30</b> 18
4.0	BAB 35 <b>40</b> 18
5.0	BAB 35 <b>50</b> 18



## Mini Ball Impression Coping

GICA
------




## Mini Ball Analog

BANL
------




## Socket Spacer

Art. No.	GBIC3L GBIC2L
----------	------------------



## Female Socket


Art. No.	BPF3 (300~500gf) BPF2 (500~700gf)
----------	--------------------------------------



※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the ball abutment with fixture.

# Magnetic Attachment [Dome Type ]

Unit: mm, Scale 1 : 1.5 / mm



MGT4520D and MKP4530D and FXS4510

## Magnetic Assay

Application	Diameter	H	Art. No.
MKP45D	Ø4.5	2.0	MGT 45 <b>20 D</b>
MKP55D	Ø5.5	2.0	MGT 55 <b>20 D</b>



## Implant Keeper Diameter Ø4.5

G/H	Art. No.
1.0	MKP 45 <b>10 D</b>
2.0	MKP 45 <b>20 D</b>
3.0	MKP 45 <b>30 D</b>
4.0	MKP 45 <b>40 D</b>
5.0	MKP 45 <b>50 D</b>
6.0	MKP 45 <b>60 D</b>



## Implant Keeper Diameter Ø5.5

G/H	Art. No.
1.0	MKP 55 <b>10 D</b>
2.0	MKP 55 <b>20 D</b>
3.0	MKP 55 <b>30 D</b>
4.0	MKP 55 <b>40 D</b>
5.0	MKP 55 <b>50 D</b>
6.0	MKP 55 <b>60 D</b>




※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the magnetic abutment with fixture.



# Magnetic Attachment [Flat Type ]

Unit: mm, Scale 1 : 1.5 / mm



MGT4515 and MKP4530 and FXS4510

## Magnetic Assay

Application	Diameter	H	Art. No.
MKP45	Ø4.5	1.5	MGT 45 <b>15</b>
MKP45	Ø4.5	2.0	MGT 45 <b>20</b>
MKP55	Ø5.5	1.5	MGT 55 <b>15</b>
MKP55	Ø5.5	2.0	MGT 55 <b>20</b>

Ø4.5  
1.5  Retention Force 400gf

Ø4.5  
2.0  Retention Force 450gf

Ø5.5  
1.5  Retention Force 700gf

Ø5.5  
2.0  Retention Force 750gf

## Implant Keeper Diameter Ø4.5

G/H	Art. No.
<b>1.0</b>	MKP 45 <b>10</b>
<b>2.0</b>	MKP 45 <b>20</b>
<b>3.0</b>	MKP 45 <b>30</b>
<b>4.0</b>	MKP 45 <b>40</b>
<b>5.0</b>	MKP 45 <b>50</b>
<b>6.0</b>	MKP 45 <b>60</b>




## Implant Keeper Diameter Ø5.5

G/H	Art. No.
<b>1.0</b>	MKP 55 <b>10</b>
<b>2.0</b>	MKP 55 <b>20</b>
<b>3.0</b>	MKP 55 <b>30</b>
<b>4.0</b>	MKP 55 <b>40</b>
<b>5.0</b>	MKP 55 <b>50</b>
<b>6.0</b>	MKP 55 <b>60</b>







※ Note: 1) It is recommended to keep the torque level at 25~30 N·cm to tighten the magnetic abutment with fixture.

# Surgical Kit [Full]




UXIF

## Kit Contents

- |                         |   |                     |                      |   |              |
|-------------------------|---|---------------------|----------------------|---|--------------|
| • First guide drill     |  | XLD 22 29           | • Parallel pin x4    |  | XPP 162220 T |
|                         |  | XLD 22 35 (Stopper) |                      |   |              |
| • Second guide drill    |  | XLD 26 29           | • Hand-piece adapter |  | XID 27 H     |
|                         |  | XLD 26 35 (Stopper) |                      |  | XID 32 H     |
| • Final drill (Stopper) |  | XFD 34 35           | • Ratchet adapter    |  | XID 26 W     |
|                         |  | XFD 38 35           |                      |  | XID 32 W     |
|                         |  | XFD 43 35           |                      |   |              |
|                         |  | XFD 48 35           |                      |   |              |
| • Final drill           |  | XFD 34 29           | • Path pin x2        |  | XMFP A2      |
|                         |  | XFD 38 29           |                      |   |              |
|                         |  | XFD 43 29           | • Hex driver         |  | XHD 26 T     |
|                         |  | XFD 48 29           |                      |   |              |
|                         |  | XFD 58 31           | • Drill extension    |  | XDE          |
| • Countersink           |  | XCS 36 29 SW        | • Depth gauge        |  | XDGL         |
|                         |  | XCS 40 29 SW        |                      |   |              |
|                         |  | XCS 45 29 SW        | • Ratchet            |  | XRCA1        |
|                         |  | XCS 50 29 SW        |                      |   |              |
|                         |  | XCS 60 29 SW        |                      |   |              |
|                         |  | XCS 70 29 SW        |                      |   |              |

# Surgical Kit [Full]



Ref. No. UXNF

Scale 0.7:1

### First Guide Drill



XLD2229



XLD2235  
(Stopper)



XLD2629



XLD2635  
(Stopper)

### Second Guide Drill



XFD 34 29



XFD 38 29



XFD 43 29



XFD 48 29



XFD 58 31

### Final Drills



XFA 27 H



XFA 32 H

### Handpiece Adapter

### Parallel Pin



XPP162220T



### Countersink Drills



XCS 3629 SW



XCS 4029 SW



XCS 4529 SW



XCS 5029 SW



XCS 6029 SW



XCS 7029SW

### Ratchet Adapter



XFA 26 W



XFA 32 W

### MF Path Pin



XMFPA 2



### Final Drills (Stopper)



XFD 34 35



XFD 38 35



XFD 43 35



XFD 48 35

### Hex Driver Drill Extension



XHD 26 T



XDE

### Depth Gauge



XDGL

### Ratchet



XRCA 1

# Surgical Kit [Stopper Drill Full]



## Stopper Drill Full Kit

UXSF

### Kit Includes

#### First Guide Drill



XLD2208S

#### Second Guide Drill



XLD2608S

#### Final Drill



XFD3408S



XFD3808S



XFD4308S



XFD4808S



XFD5831

#### Hand-piece Adapter



XFA27H



XFA32H



XLD2210S



XLD2610S



XFD3410S



XFD3810S



XFD4310S



XFD4810S

#### Ratchet Adapter



XFA26W



XFA32W



XLD2212S



XLD2612S



XFD3412S



XFC3812S



XFD4312S



XFD4812S

#### HexDriver Drill Extension



XHD26T



XDE

#### Parallel Pin



XPP162220T



XPP162220T



XFC3607



XFC4007



XFC4507



XFC5007



XFC6007



XFC7007

#### Countersink

#### Hand-piece Ratchet



XHD25H



XHD25W

#### MF Path Pin



XMFPA2

#### Depth Gauge



XDGL

#### Torque Wrench




XNTW




# Digital Guide | Full Kit for FX 36/40/45

Unit: mm, Scale 1 : 1



XGSFK

## SuperLine FSX4508 Guide Drilling Sequence



























# Surgical Kit [Standard ]




UXIFN

## Kit Contents
















• First guide drill		XLD 22 31	• Parallel pin		XPP 162220 T
• Second guide drill		XLD 26 31			XPP 162226 T
• Final drill		XFD 34 31	• Hand-piece adapter		XID 30 H
		XFD 38 31	• Ratchet adapter		XID 26 W
		XFD 43 31	• MF path pin x2		XMFPA2
		XFD 48 31	• Hex driver		XHD 26 T
		XFD 58 31			XHD 25 H
• Countersink		XCS 36 29 SW	• Drill extension		XDE
		XCS 40 29 SW	• Ratchet		XRCA1
		XCS 45 29 SW			
		XCS 50 29 SW			
		XCS 60 29 SW			
		XCS 70 29 SW			

# Surgical Kit [Short Implant ]




XSIK

## Kit Contents

• First guide drill		XLD <b>22</b> 07	• Parallel pin × 2		XPP 162212
• Second guide drill		XLD <b>26</b> 07	• Hand-piece adapter		XID <b>27</b> H
• Final drill		XFC <b>36</b> 07	• Ratchet adapter		XID <b>26</b> W
		XFC <b>40</b> 07	• MF path pin x2		XMFP A2
		XFC <b>45</b> 07	• Hex driver		XHD <b>26</b> T
		XFC <b>50</b> 07			XHD <b>25</b> H
		XFC <b>60</b> 07	• Ratchet		XRCA1
		XFC <b>70</b> 07			

# Drill Stopper Kit



XDS

## Kit Contents

- Guide drill stopper / First, Second

 XLD 22 35

 XLD 26 35

- Final drill stopper / 34, 38

 XFD 34 35

 XFD 38 35

- Final drill stopper / 43, 48

 XFD 43 35

 XFD 48 35

- Stopper-first guide drill, second guide drill

XLDST 14   XLDST 12   XLDST 10   XLDST 08   XLDST 06



- Stopper- final drill / 34, 38

XFDST 14   XFDST 12   XFDST 10   XFDST 08   XFDST 06



- Stopper- final drill / 43, 48

XFDST 14L   XFDST 12L   XFDST 10L   XFDST 08L   XFDST 06L




# Drill



Unit: mm, Scale 1 : 1 / mm


## Guide Drill

Diameter	L	Art. No.
Ø1.5	20	XGD 15 <b>20</b>
Ø1.5	28	XGD 15 <b>28</b>




## First Drill

Diameter	L	Art. No.
Ø2.0	29	XFD 20 <b>29</b>
Ø2.0	33	XFD 20 <b>33</b>




## First Guide Drill

Diameter	L	Art. No.
Ø2.2	26	XLD 22 <b>07</b>
Ø2.2	29	XLD 22 <b>29</b>
Ø2.2	31	XLD 22 <b>31</b>
Ø2.2	33	XLD 22 <b>33</b>
Ø2.2	35	XLD 22 <b>35</b>
Ø2.2	41	XLD 22 <b>41</b>



## Second Guide Drill

Diameter	L	Art. No.
Ø2.6	26	XLD 26 <b>07</b>
Ø2.6	29	XLD 26 <b>29</b>
Ø2.6	31	XLD 26 <b>31</b>
Ø2.6	33	XLD 26 <b>33</b>
Ø2.6	35	XLD 26 <b>35</b>
Ø2.6	41	XLD 26 <b>41</b>



※ Note: Drill speed 1,000rpm, 30~45N·cm with irrigation.

# Drill

Unit: mm, Scale 1 : 1 / mm


## Pilot Drill

Diameter	L	Art. No.
Ø3.0	30	XPD 20 30




## Final Drill | For Short Implant

Diameter	L	Art. No.
Ø3.6	27	XFC 36 07
Ø4.0	27	XFC 40 07
Ø4.5	27	XFC 45 07
Ø5.0	27	XFC 50 07
Ø6.0	27	XFC 60 07
Ø7.0	27	XFC 70 07




## Final Drill | Length-29mm

Diameter	L	Art. No.
Ø2.9	29	XFD 34 29
Ø3.35	29	XFD 38 29
Ø3.85	29	XFD 43 29
Ø4.4	29	XFD 48 29
Ø5.4	29	XFD 58 29SW



## Final Drill | Length-31mm

Diameter	L	Art. No.
Ø2.9	31	XFD 34 31
Ø3.35	31	XFD 38 31
Ø3.85	31	XFD 43 31
Ø4.4	31	XFD 48 31
Ø5.4	31	XFD 58 31



※ Note: Drill speed 1,000rpm, 30~45N·cm with irrigation.




# Drill

Unit: mm, Scale 1 : 1 / mm


## Final Drill | Length-33mm

Diameter	L	Art. No.
Ø2.9	33	XFD 34 33
Ø3.35	33	XFD 38 33
Ø3.85	33	XFD 43 33
Ø4.4	33	XFD 48 33




## Final Drill | Length-35mm | Stopper

Diameter	L	Art. No.
Ø2.9	35	XFD 34 35
Ø3.35	35	XFD 38 35
Ø3.85	35	XFD 43 35
Ø4.4	35	XFD 48 35



## Final Drill | Length-41mm

Diameter	L	Art. No.
Ø2.9	41	XFD 34 41
Ø3.35	41	XFD 38 41
Ø3.85	41	XFD 43 41
Ø4.4	41	XFD 48 41




※ Note: Drill speed 1,000rpm, 30~45N·cm with irrigation.

# Drill

Unit: mm, Scale 1 : 1 / mm


## Harvest Drill | Length-35mmStopper

Diameter	L	Art. No.
Ø2.85	35	XFH 34 35
Ø3.3	35	XFH 38 35
Ø3.85	35	XFH 43 35
Ø4.4	35	XFH 48 35



## Countersink

Diameter	L	Art. No.
Ø3.6	29	XCS 36 29 SW
Ø4.0	29	XCS 40 29 SW
Ø4.5	29	XCS 45 29 SW
Ø5.0	29	XCS 50 29 SW
Ø6.0	27	XCS 60 29 SW
Ø7.0	27	XCS 70 29 SW



## Round Bur

Diameter	L	Art. No.
Ø2.0	33	XRB 20 33
Ø3.0	33	XRB 30 33



※ Note: Drill speed 1,000rpm, 30~45N·cm with irrigation.

# Stopper

Unit: mm, Scale 1 : 1 / mm

## Stopper | For first guide drill, second guide drill

Drilling Depth	L	Art. No.
14	4.6	XLDST 14
12	6.6	XLDST 12
10	8.6	XLDST 10
08	10.6	XLDST 08
06	12.6	XLDST 06
04	14.6	XLDST 04
02	16.6	XLDST 02



## Stopper | For final drill 3435, 3835

Drilling Depth	L	Art. No.
14	4.6	XFDST 14
12	6.6	XFDST 12
10	8.6	XFDST 10
08	10.6	XFDST 08
06	12.6	XFDST 06
04	14.6	XFDST 04
02	16.6	XFDST 02



## Stopper | For final drill 4335, 4835

Drilling Depth	L	Art. No.
14	4.6	XFDST 14L
12	6.6	XFDST 12L
10	8.6	XFDST 10L
08	10.6	XFDST 08L
06	12.6	XFDST 06L
04	14.6	XFDST 04L
02	16.6	XFDST 02L



# Instrument


Unit: mm    Scale 1:1

## Fixture Adapter | Hex 2.5mm


Type	Length	Art. No.
Handpiece	27	XFA 27 H
	30	XFA 30 H
	32	XFA 32 H
Ratchet	21	XFA 21 W
	24	XFA 24 W
	26	XFA 26 W
	29	XFA 29 W
	32	XFA 32 W

- Use the handpiece adapter to transfer the implant fixture.
- **NOTE:** Ratchet Adapters are intended for final adjustments of the fixtures after placement of the fixture by the handpiece adapter. Ratchet adapters are not intended for transferring the fixture to or from the osteotomy site.
- To reduce the risk of accidental swallowing or aspiration, use with caution when transferring fixtures and its related restorative components

### Handpiece Adapter




### Ratchet Adapter




## Parallel Pin | For first guide drill, second guide drill

Diameter	Art. No.
Ø1.6	XPP 162212
Ø1.6	XPP 162220T
Ø1.6	XPP 162226T



## Path Pin

L	Art. No.
18.6	XMFP A2






# Instrument

Unit: mm, Scale 1 : 1 / mm


## Hex Driver | Hex 1.28mm

Type	L	Art. No.
Hand-piece	<b>25</b>	XHD <b>25</b> H
Ratchet	<b>21</b>	XHD <b>21</b> W
	<b>25</b>	XHD <b>25</b> W
	<b>27</b>	XHD <b>27</b> W
Manual	<b>26</b>	XHD <b>26</b> T
	<b>30</b>	XHD <b>30</b> T




## Slot Driver

Type	Art. No.
Ratchet	SDA <b>17</b> R
	SDA <b>25</b> R




## Drill Extension

XDE
-----



## Driver | Manual

Type	Art. No.
Manual	XHDHT




# Instrument

Unit: mm

**Adapter** | For screw & ball abutment | Scale 1 : 1 / mm

Type	Art. No.
Hand-piece	XMAA1
Ratchet	XMA 21W



**Adapter** | For mini ball abutment | Scale 1 : 1 / mm

IPST21W
---------



**Ratchet**

XRCA1
-------



**Torque Wrench** | Scale 1 : 0.7 / mm

XNTW
------



**Depth Gauge**

XDGL
------

※ Note: One side of Depth Gauge measures the osteotomy depth and the other side measures the gingival height from the top of the implant.



**Tissue Punch** | Scale 1 : 1 / mm

XTS40
-------



※Punching size : Ø4.0



# DASK [Dentium Advanced Sinus Kit]



## Kit Contents

- DASK drill



- Stopper



- Sinus elevation instrument



# Sinus Bur Kit / Sinus Kit



SDK

## Kit Contents

- DASK drill



XRT**33**2035



XRT**06**4025



XRT**37**2035



XRT**08**4025



XED**33**1035D



XST**08**3025

- Stopper

XFDST **08**



XFDST **06**




XFDST **04**



XFDST **02**



## Sinus Kit



XSKL


# DASK / Sinus Bur Kit / Sinus Kit

Unit: mm

DASK Drill | Scale 1 : 1.2 / mm

Type	DASK Drill #	Art.No.
Crestal Approach	DASK Drill #1	XRT <b>33</b> 2035
	DASK Drill #2	XRT <b>37</b> 2035
	DASK Drill #3	XED <b>33</b> 1035D
Lateral Approach	DASK Drill #4	XRT <b>06</b> 4025
	DASK Drill #5	XRT <b>08</b> 4025
	DASK Drill #6	XST <b>08</b> 3025

※ Note: Drill speed 800 to 1,200rpm, 30~45N·cm with irrigation.




Stopper | For XRT332035, XRT372035, XED331035D | Scale 1 : 1 / mm

Drilling Depth	L	Art.No.
<b>08</b>	10.6	XFDST <b>08</b>
<b>06</b>	12.6	XFDST <b>06</b>
<b>04</b>	14.6	XFDST <b>04</b>
<b>02</b>	16.6	XFDST <b>02</b>



Sinus Elevation Instrument | Scale 1 : 0.45 / mm




# DASK Simple




## Crestal approach

### Bone chip lifting & Compaction

Safe, easy and efficient surgery possible  
without perforation



50rpm 70N.cm  
Non irrigation




Sinus Lifting




Sinus Lifting




Round Instrument



Bone Spreader



# DETAILED SPECIFICATION



## Crestal Drill



Top view




Diameter  
Ø 3.0

※ Recommended usage frequency 50 times

Diameter	L	Art. No
Ø 3.0	35.0mm	SRT3035AS
Ø 3.5	35.0mm	SRT3535AS


## Lateral Drill



※ Recommended usage frequency 40 times

Diameter	L	Art. No
Ø 4.0	28.0mm	SLD0428
Ø 4.0	40.0mm	SLD0440

## Bone Spreader



Diameter	L	Art. No
Ø 5.0	33.5mm	SSP33


# DASK Simple

## Lateral approach




✓ Improvement in the user-friendliness and easy-operation

Lateral Drill with the ability to freely shape the Sinus Lateral from various angles




✓ Easy access for both vertical and lateral approach

1,000rpm 30~45N.cm with irrigation



DASK  
Ø6.0~8.0




DASK simple  
Ø4.0





# Ridge Spreader Kit



## XRSK

### Bone Chisel

XBC305013
-----------



### Ratchet

XRCA1
-------



[Unit: mm, Scale 0.6 : 1]

### Ridge Spreader Drills

Diameter	L	Art No.
Ø1.4 / Ø2.4	35	RS142435
Ø2.0 / Ø3.2	35	RS203235
Ø2.6 / Ø3.6	35	RS263635



### Round Bur

Diameter	L	Art No.
Ø4.0	35	XR4035



### Ratchet Adapter

XRA3917
---------



### Mini Saw

Diameter	L	Art No.
Ø8.0	25	XDS8025



[Unit: mm, Scale 1 : 1]

# Osteotome Kit


Unit: mm

## Osteotome Kit

Osteotome compresses the bone laterally, providing denser bony interface rather than removing valuable bone from the surgical site.




### Osteotome | Final drill type | Scale 1 : 0.4 / mm



## Osteotome Kit

Type	Art.No.	Ø A	Ø B
<b>XOFK</b> Type <b>A</b> (Convex)	XOF 20 <b>A1</b>	Ø 1.7	Ø 2.8
	XOF 34 <b>A1</b>	Ø 2.3	Ø 2.8
	XOF 38 <b>A1</b>	Ø 2.7	Ø 3.2
	XOF 43 <b>A1</b>	Ø 2.8	Ø 3.8
	XOF 48 <b>A1</b>	Ø 3.0	Ø 4.3
<b>XOFBK</b> Type <b>B</b> (Concave)	XOF 20 <b>B1</b>	Ø 1.7	Ø 2.8
	XOF 34 <b>B1</b>	Ø 2.3	Ø 2.8
	XOF 38 <b>B1</b>	Ø 2.7	Ø 3.2
	XOF 43 <b>B1</b>	Ø 2.8	Ø 3.8
	XOF 48 <b>B1</b>	Ø 3.0	Ø 4.3



XOFK (Type A)


XOFBK (Type B)

# Trephine Kit

## Trephine Bur


- Excellent fine cutting
- Strong engagement when attaching the trephine to cortical bone
- Cut-outs facilitates ease of harvest retrieval
- 5 scale marks on the Trephine drill from 2mm to 10mm
- Easy harvesting

Unit: mm




**Trephine Bur** | Scale 1 : 0.5 / mm

Outer Diameter	Inner Diameter	Art. No.
Ø <b>3.0</b>	Ø <b>2.3</b>	XTP 24 03
Ø <b>4.0</b>	Ø <b>3.3</b>	XTP 34 04
Ø <b>5.0</b>	Ø <b>4.2</b>	XTP 44 05
Ø <b>6.0</b>	Ø <b>5.2</b>	XTP 54 06
Ø <b>7.0</b>	Ø <b>6.1</b>	XTP 64 07
Ø <b>8.0</b>	Ø <b>7.1</b>	XTP 74 08
Ø <b>9.0</b>	Ø <b>8.0</b>	XTP 84 09
Ø <b>10.0</b>	Ø <b>9.0</b>	XTP 94 10




## Trephine Kit

XIT



# Prosthetic Kit



XIP

**Hex Driver S/T** | Scale 1 : 1 / mm

XHD 15



**Hex Driver L/T** | Scale 1 : 1 / mm



**Slot Driver** | Scale 1 : 1 / mm

SDA 25 R



**Mount Adapter** | Scale 1 : 1 / mm

XMAA1



**Hex Driver T/W** | Scale 1 : 1 / mm

XHD 25 W



**T/W Adapter** | Scale 1 : 1 / mm

XMA 21 W



**Mini Ball Adapter** | Scale 1 : 1 / mm

IPST 21 W



**Torque Wrench** | Scale 1 : 0.5 / mm

XNTW



# Planning Kit

Unit: mm, Scale 1 : 1 / mm



## Kit Contents

### Diameter Ø4.5 | Combi & Dual abutment

G/H	Art.No.
1.5	PDAB 45 <b>15</b>
2.5	PDAB 45 <b>25</b>
3.5	PDAB 45 <b>35</b>
4.5	PDAB 45 <b>45</b>
5.5	PDAB 45 <b>55</b>



### Diameter Ø5.5 | Combi & Dual abutment

G/H	Art.No.
1.5	PDAB 55 <b>15</b>
2.5	PDAB 55 <b>25</b>
3.5	PDAB 55 <b>35</b>
4.5	PDAB 55 <b>45</b>
5.5	PDAB 55 <b>55</b>



### Diameter Ø6.5 | Combi & Dual abutment

G/H	Art.No.
1.5	PDAB 65 <b>15</b>
2.5	PDAB 65 <b>25</b>
3.5	PDAB 65 <b>35</b>
4.5	PDAB 65 <b>45</b>
5.5	PDAB 65 <b>55</b>



### Angled 15° | Angled abutment

Diameter	G/H	Art.No.
Ø4.5	2.0	PAAB 15 45 <b>20</b>
Ø4.5	4.0	PAAB 15 45 <b>40</b>
Ø5.5	2.0	PAAB 15 55 <b>20</b>
Ø5.5	4.0	PAAB 15 55 <b>40</b>



### Angled 25° | Angled abutment

Diameter	G/H	Art.No.
Ø4.5	2.0	PAAB 25 45 <b>20</b>
Ø4.5	4.0	PAAB 25 45 <b>40</b>
Ø5.5	2.0	PAAB 25 55 <b>20</b>
Ø5.5	4.0	PAAB 25 55 <b>40</b>




# Prosthetic and Laboratory Instrument

Unit: mm

Hex Driver | Hex 1.28 mm | Scale 1 : 1 / mm

L	Art. No.
13	XHD <b>13</b>
15	XHD <b>15</b>
21	XHD <b>21</b>
28	XHD <b>28</b>



Reamer Guide for Combi/Dual Abutment | Scale 1 : 1 / mm

Diameter	Art. No.
Ø <b>4.5</b>	CRG <b>45</b> L
Ø <b>5.5</b>	CRG <b>55</b> L
Ø <b>6.5</b>	CRG <b>65</b> L



Reamer Guide for Screw Abutment | Scale 1 : 1 / mm

Type	Art. No.
<b>Bridge</b>	SRG <b>BL</b>
<b>Single</b>	SRG <b>SL</b>



# Prosthetic and Laboratory Instrument

Reamer Handle | Scale 1 : 0.5 / mm

CRH



Hand Wrench | Scale 1 : 1 / mm

XHW



Reamer (Combi/Dual Abutment) | Scale 1 : 1 / mm

CRM



Reamer (Screw Abutment) | Scale 1 : 1 / mm


SRM





# Help Kit

Unit: mm



XIH


## Cover & Abutment Screw Remover | Scale 1 : 1 / mm

L	Art. No
25	XRRHD




## Abutment Hex Remover | Scale 1 : 1 / mm

L	Art. No
20	XRHR20
25	XRHR25





## Screw Remover | Scale 1 : 1 / mm

L	Art. No
25	XRFS2S
33	XRFS2



## Fixture Remover | Scale 1 : 1 / mm

Type	Art. No
Remover	XRFRT
	XRFRTF
	XRFRTTO
Wrench	XRFRW



## Screw Tap Repair | Scale 1 : 1 / mm

Type	Art. No
Tap	XRSTR
11° Guide	XRSG11
8° Guide	XRSG8






# SURGICAL **MANUAL**


Surgical Drill Sequence I	74
Surgical Drill Sequence II	76
Drilling Depth Guide	78
Fixture Connection	80
Installation Procedure & Warnings	81
Surgical Kit Maintenance	82



# Surgical Drill Sequence I




## Drilling Sequence Guide (Final Drill)



**During Fixture Insertion, 30 ~ 45N·cm Torque at 20rpm is Recommended**

- Countersink drill is used in cases with dense cortical bone.
- If the bone density is D1~D3, it is recommended to countersink after final drill.
- The actual diameter of the Countersink drill is 0.1mm larger than the fixture platform.

**Determination of Fixture Top Level**



Top level of fixture needs to be located 0.5mm below the marginal crestal bone level to minimize bone loss after implantation. However, only for the fixture of 7mm length, top level of fixture should be located 0.5mm above the marginal crestal bone level.



**Depth Indication**



- Use the depth gauge after first drill / First guide drill to check depth of drilling.
- Place the depth gauge against the wall of the osteotomy.

# Surgical Drill Sequence II

## Drilling Sequence Guide (Havrest Drill)




△ : During the 4.3/4.8 fixture insertion into the bone density of D3~D4, the 3.4/3.8 harvest drilling process can be skipped.


Harvest Drill

Simultaneous and effective autogenous bone collection during the final drilling procedure using a specially designed drill


- Sharp pointed tip design of the drill prevents drill chatter and helps guide drill path.
- Available drill stopper helps control drill depth for safe and efficient bone collection, especially in the buccal side of ridge.
- Recommended drill speed of less than 100rpm/50N-cm helps preserve the vital autogenous bone.
- Excellent clinical results may be achieved when harvested autogenous bone is combined with OSTEON™ II.



Bone Collection in the Buccal Side of Ridge: 50~200rpm / 30~50N-cm



# Drilling Depth Guide



## Countersink Depth Guide

- Drilling depth of the countersink depends on the patient's bone quality.
- If the bone density is D1~D2, it is recommended to drill up to the top line ( I ) of laser mark on the countersink.
- If the bone density is D3~D4, it is recommended to drill up to the bottom line ( II ) of laser mark on the countersink.

Platform: Ø3.6 / Body: Ø3.4 (1000rpm / 30~45N·cm)

First guide drill

Second guide drill

Final drill Ø3.6

Countersink Ø3.7

Platform: Ø4.0 / Body: Ø3.8 (1000rpm / 30~45N·cm)

First guide drill

Second guide drill

Final drill Ø3.6

Final drill Ø4.0

Countersink Ø4.0

Platform: Ø4.5 / Body: Ø4.3 (1000rpm / 30~45N·cm)

First guide drill

Second guide drill

Final drill Ø3.6

Final drill Ø4.0

Final drill Ø4.5

Countersink Ø4.5





# Fixture Connection



Caution\_ When opening the fixture pack, hold the fixture container upward and engage the adapter into the fixture.




By hand-piece  
20rpm / 35N-cm



By ratchet



## Directions Using the Hand-piece / Ratchet Adapter




Hand-piece  
Adapter

Ratchet  
Adapter

The Hand-piece Adapter/Ratchet Adapter must be connected firmly together with the internal hex inside the fixture

# Installation Procedure & Warnings

## Cover Screw



By Hex Driver



Cover Screw (CS36)  
connection

## Healing Abutment



By Hex Driver



Healing Abutment connection



Healing Abutment (HAB402020L)  
connection in thin gingiva

## Warnings

Dental Implant surgery and restoration involve complex dental procedures. Appropriate and adequate training in proper technique is strongly recommended prior to use.

- Improper medical examination and/or treatment plan can result in implant failure and/or loss of supportive bone.
- Improper initial stability and/or excessive occlusal forces during healing period may lead to osseointegration failure.
- Excessive insertion torque may lead to mechanical failure or implant biologic failure due to bone compression and necrosis.
- When forces or loads are greater than its design, implant or abutment fracture could happen. Therefore clinicians should make careful decisions with regards to clinical treatment planning to minimize the risk of fracture. Appropriate implant quantity, occlusal interface and a nightguard are essential. Potential excessive loading conditions may include the following:

- 01 Inadequate number of implants are placed.
- 02 Implant width and/or length are inappropriate for a treatment site.
- 03 Prosthesis which has excessive cantilever length due to inadequate biomechanical design
- 04 Continuous occlusal force are generated by incomplete connection between implant and abutment and/or abutment screw loosening.
- 05 Direct Casting Abutment angles are greater than 30° from the vertical axis of the implant.  
Direct Abutments are not for angulation.
- 06 Occlusal interferences causing excessive lateral forces
- 07 Patient parafunctions such as bruxism
- 08 Inadequate dental laboratory casting procedures
- 09 Improper prosthesis fit
- 10 Trauma from patient habits or accidents
- 11 Excessive marginal bone loss caused by inadequate bone width and/or advanced periimplantitis

# Surgical Kit Maintenance

## Manual Cleaning and Sterilization Procedure

It is important to use protective clothing and face shield while cleaning contaminated instruments. Always wear protective glasses, mask, gloves, etc. for your safety.

### Cleaning

- 1 Rinse instruments immediately after use under running tap water (<40°C) for a minimum of one (1) minute to remove all debris including extraneous body fluids, bone debris and tissue.
- 2 Soak all instruments immediately after rinsing in an enzymatic cleaning solution\* for 10 to 20 minutes (Do not soak overnight).  
  
\* Follow manufacturer’s instructions and observe recommended cleaning solution concentrations (enzymatic detergent with a pH level between 7-10 and temperature not to exceed 40°C). Do not use incompatible cleaning solutions to clean instruments.
- 3 For internal irrigation drills, use a 1mL syringe and a 25 gauge needle to clean the drill irrigation hole with a minimum of 0.2 mL of the prepared cleaning solution. Repeat this step two (2) more times for a total of three (3) rinses.
- 4 Scrub with a soft brush for a minimum of 1 (one) minute to remove any debris inside the drill irrigation hole.
- 5 Rinse the instruments under running tap water (<40°C) for a minimum of 1 minute. Use a 1mL syringe and a 25 gauge needle with a minimum of 0.2 mL of tap water to forcefully flush inside the drill irrigation hole. Repeat flushing of drill irrigation hole two (2) more times for a total of three (3) flushings.
- 6 Place instruments into an ultrasonic cleaner with neutral detergent\*\*. Keep instruments inside the ultrasonic bath for 15 minutes using a frequency of 25-50 kHz. Ensure multiple instruments placed within the bath remain separated.  
  
\*\* Follow manufacturer’s instructions and observe recommended neutral detergent solution concentrations (neutral detergent with a pH level between 7-10 and temperature not to exceed 40°C). Do not use incompatible neutral detergent solutions to clean instruments.
- 7 Rinse instruments thoroughly with running tap water (<40°C) for a minimum of 1 (one) minute until all traces of neutral detergent solution are removed. Rinse inside drill irrigation hole using a 1mL syringe and a 25 gauge needle with a minimum of 0.2 mL of tap water. Repeat rinsing drill irrigation hole two (2) more times for a total of three (3) rinses.
- 8 Gently wipe instruments with a soft lint-free cloth or place the instruments in a drying cabinet (60°C for less than 10 hours) until fully dry. Blow residual water from drill irrigation hole using a 1mL syringe and a 25 gauge needle. Visually inspect instruments in a well-lit area to ensure they are clean, dry and free of residue.
- 9 Clean instrument trays with a germicidal cleaner prior to returning instruments into Kit.
- 10 Always check for damage or corrosion after rinsing and drying.

### Sterilization

Dentium recommends either the Pre-vacuum or Gravity autoclave methods for sterilization under the conditions described below. However, autoclave performance can affect the efficacy of this process. Healthcare facilities should validate their sterilization processes employing the actual equipment and operators that routinely sterilize instruments.

All autoclaves/sterilizers should be regularly validated, maintained and checked in accordance with EN 285/EN 13060, EN ISO 17665, ANSI AAMI ST79 to ensure compliance with these and related standards. Make sure packaging is suitable for steam sterilization.

#### Recommended Sterilization Parameters

Method-Moist Heat Sterilization	Pre-vacuum	Gravity
Set Point Temperature	132 °C	132 °C
Exposure time	4 minutes	30 minutes
Drying time	20 minutes	40 minutes



# PROSTHESIS **MANUAL**

## **Prosthetic Introduction**

Understanding the Implant and Prosthesis	84
Types of Abutment	85
Dual Abutment	86
Combi Abutment	87
Dual Milling / Anged / Temporay / Metal-Casting Abutment	88
Screw Abutment / Angled Screw Abutment	90
Points to Consider in Abutment Selection	91
Minimum Height Requirement for SuperLine Prosthetic Abutment	92

## **Prosthetic Procedure 1**

Abutment Level- Dual Abutment	94
Abutment Level- Combi Abutment	97

## **Prosthetic Procedure 2**

Fixture Level [Pick-up Type]- Dual Abutment	100
Fixture Level [Transfer Type]- Dual Abutment	103
Fixture Level [Transfer Type]- Dual Milling Abutment	106
Fixture Level [Pick-up Type]- Angled Abutment	108
Fixture Level- Metal-Casting Abutment	111
Fixture Level [Pick-up Type]- Temporary Abutment	112

## **Prosthetic Procedure 3**

Abutment Level [Transfer Type]- Screw Abutment	113
	114

## **Cementation Repair Method (SCRIP)**


117

## **Prosthetic Procedure 4**

Positioner	119
Ball Attachment	120
Magnetic Attachment	122
	123



# Understanding the Implant and Prosthesis




## Biological Connection

- The tapered conical hex connection between implant and abutment interface provides hermetic sealing.
- The biological connection distributes the load to the fixture evenly. Therefore it may minimize bone loss.
- All implant diameters share the same internal connection. One abutment screw fits all abutments and fixtures.







## Types of Abutment (Abutments are available in various diameters & gingival heights)

- Dual Abutment
- Combi Abutment
- Dual Abutment
- Dual Milling Abutment
- Angled Abutment (15°/25°)
- Metal-Casting Abutment
- Temporary Abutment (Plastic & Titanium)
- Screw Abutment
- Angled Screw Abutment (15°/ 30°)
- Positioner Attachment
- Ball Attachment
- Magnetic Attachment





# Types of Abutment

One-Piece	Two-Pieces	
 Combi Abutment	 Hex	 Non-hex
	 Hex	 Non-hex
	 Hex	 Non-hex
	 Hex	 Non-hex
 Screw Abutment	 Hex	 Non-hex
	 Hex	 Non-hex
	 Hex	 Non-hex
 Cylinder	 Hex	 Non-hex
Abutment Level	Fixture Level	

- Straight abutments are Dual and Combi Abutment.
- Depending on the insertion angle and position of the fixture, the Angled or / Metal - Casting Abutment may be used.
- The Screw Abutment can be used when prosthesis retrieval is anticipated.

## Selection Guideline

Ideal emergence profile for each tooth



# Dual Abutment



- It is possible to take an impression at both fixture level and abutment level.  
(A Dual Abutment may be interchanged with a Combi Abutment)
- For abutment level impressions, the same prosthetic procedures apply to both Dual and Combi Abutments.
- For fixture level impressions, the abutment selection takes place on the master model.
- For fixture level impressions, a precise positioning jig for abutment may be required.
- Either hex or non-hex abutments may be used, according to operator's preference.

\* If a cement retained restoration requires retrieval, cutting a hole in the occlusal surface would allow access to the screw to permit removal.

## Hex / Non-hex


	Hex	Non-hex
Positioning Jig	Unnecessary	Required
Radiograph	Required	Unnecessary

## Dual Abutment (Hex / Non-hex)

Diameter	G/H	Vertical Angle
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	5°
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	6°
Ø6.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	7°



# Combi Abutment



- The Combi Abutment is used when the implant position is optimal.
  - If the abutment selection is made in the mouth, gauge the thickness of mucosa with the depth gauge to measure the gingival height thus allowing the appropriate abutment height.
  - The Impression is taken with the snap cap.
  - When using the Combi Abutment, it remains in the mouth after the impression is taken. (Do not remove or change its position)
  - Tighten abutment screw to 25 - 35 N·cm. (retighten again before seating final prosthesis)
- \* If the Combi Abutment is too long it can be adjusted 1.5mm to the bottom of the laser mark on the vertical stack of the abutment. The Combi Abutment has a short analog for the 1.5mm adjustment.
- \* A resin jig can be made to record the reduction if reduced more the 1.5mm.

## Combi Abutment Line Up

Diameter	G/H	Vertical Angle
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	5°
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	6°
Ø6.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	7°



# Dual Milling / Angled / Temporary / Metal-Casting Abutment




## Dual Milling Abutment

- Impression is taken at fixture level.
- When using a non-hex abutment a precise seating jig should be used.
- Either hex or non-hex abutments may be used, according to operators preference.

\* If a cement retained restoration requires retrieval, cutting a hole in the occlusal surface would allow access to the screw for removal.

## Angled Abutment

- The Angled Abutment is recommended when the restoration path of insertion is unfavorable in either anterior or posterior sites.
- Retention force can be increased through milling process.



# Dual Milling / Angled / Temporary / / Metal-Casting Abutment

## Metal-Casting Abutment

- Equivalent results for a fraction of the price
- Our highly affordable metal alloy replaces expensive gold to alleviate financial burden to all.

## Temporary Abutment

- Temporary Abutments are available with titanium or plastic.
- The titanium abutment comes in both hex and non-hex with a gingival height of 1.0mm.
- The plastic abutment comes in diameters (Ø4.5, 5.5, 6.5) with a gingival height of 2.0mm.

Fixture Level Abutment (Hex / Non-hex)			
Abutment	Diameter	G/H	Angle
Dual Milling 	Ø4.0	1.0mm	X
	Ø4.5	1.5mm	
	Ø5.5	1.5 / 2.5mm	
	Ø6.5	1.5 / 2.5 / 3.5mm	
	Ø7.5	2.5 / 3.5mm	
Angled 	Ø4.5	1.5mm 2.5mm 3.5mm	15° / 25°
	Ø5.5	1.5mm 2.5mm 3.5mm	15° / 25°
Metal-Casting Abutment 	Ø4.5	1.0mm	X
Ti-Temporary 	Ø4.5	1.0mm	X
Plastic Temporary 	Ø4.5 Ø5.5 Ø6.5	2.0mm	X

# Screw Abutment



Screw Abutment



Angled Screw Abutment

**If prosthesis repair is anticipated, use of a Screw Abutment retained prosthesis enables easy retrieval.**

- Useful for connecting multiple units or when there is a preference for a screw retained prosthesis.
- Useful when respective long axes of implants differ. Each side tapers by 30° and this permits up to 60° divergence between two abutments.
- Useful when the prognosis of an adjacent restoration is not ideal thus permitting easy retrieval and modification of the restoration.

## Ti-Retaining Screw (1.8mm - body diameter)

- Can minimize screw loosening due to increased approximal space.
- Can endure various kinds of masticatory force.



## Screw Abutment

Diameter	G/H
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm

## Angled Screw Abutment

Diameter	G/H	Angled
Ø4.5	1.0mm	15°
Ø5.5	1.5mm	30°





# Points to Consider in Abutment Selection

## Considerations in Selecting an Abutment

- Esthetic requirement
- Implant angulation
- Implant location
- Fixture installation depth (Gingival height)
- Interarch distance
- Prosthesis type
- Dentist & dental technician's preference

## Impression of Implant


According to the case the impression can be taken at abutment or fixture level.

Fixture Level	Abutment Level
1. Dual Abutment	1. Dual Abutment
2. Dual Milling Abutment	2. Combi Abutment
3. Angled Abutment (15° / 25°)	3. Screw Abutment
4. Metal-Casting Abutment	4. Angled Screw Abutment (15° / 30°)
5. Temporary Abutment (Plastic & Titanium)	

## Abutment Impression Recommendation

Dual Abutment	Cementation type, screw-cementation type	Fixture level impression or abutment level impression
Combi Abutment	Cementation type	Abutment level impression
Angled Abutment	Cementation type, screw-cementation type	Fixture level impression
Screw Abutment	Screw retained type	Abutment level impression
Metal-Casting Abutment	Cementation type, screw-cementation type	Fixture level impression
Dual Milling Abutment	Cementation type, screw-cementation type	Fixture level impression

# Minimum Height Requirement for SuperLine Prosthetic Abutment



※ Diagram above indicates the minimum height required for SuperLine prosthetic abutment.

## Maximum Amount of Reduction Allotted for SuperLine

### Combi Abutment

- Eliminate 3.0mm from the top level Combi Abutment (laser marking:1.5mm)  
Caution \_ Damage may be caused to the screw if the abutment is reduced to less than 2.5mm above the gingival height.

### Dual Abutment

- Preparation of the abutment top is possible as follows.

Gingival Height	Preparable Amount
1.5mm	2.0
2.5mm	3.0
3.5mm	4.0
4.5mm	5.0
5.5mm	6.0

### Angled Abutment & Dual Milling Abutment

- Required minimum abutment height: at least 5.0mm above the Fixture top

### Metal-Casting Abutment

- Required minimum abutment height: at least 5.5mm above the Fixture top.


### Screw Abutment

- The Screw Abutment cannot be modified, however the Casting Abutment can be modified for interarch distance, taking reduction into consideration of the height of the retaining screw.

# Prosthetic Procedure 1

Impression Technique and Restoration Selection

## Dual / Combi Abutment




# Abutment Level- Dual Abutment

[Multiple Units]

## Clinical Procedure




## Chairside




# Abutment Level- Dual Abutment

[Multiple Units]

## Laboratory Procedure



## LabSide



# Abutment Level- Dual Abutment

[Multiple Units]



Trimming of the extended margin by using the rubber wheel



Metal framework and reamer



Reamer is used to eliminate "Lip" caused by 'snap-on' mechanism.



Metal framework after removal of "Lip"



Metal framework



Porcelain build-up

**SCRIP:** Once an access hole has been created, it can be converted to a SCRIP (Screw & Cemented Retained Prosthesis).



Final prosthesis




Access hole is made when burn-out cylinder is used to do the wax-up.



Extended margin around the metal framework due to 'snap-on' mechanism



Trim extended margin by rubber wheel



Metal framework and reamer



Eliminate the lip remnant caused by 'snap-on' mechanism by reamer.



Metal framework after removal of "Lip"



Metal framework



Final prosthesis



# Abutment Level- Combi Abutment

[Multiple Units]

## Chairside



Second stage surgery (uncovering)



Following the 2nd stage surgery, soft tissue is healed around the Healing Abutment. Healing Abutment should be selected according to the size of abutment.



Choose abutment with gingival height then tighten it to 25~30N-cm. Re tighten after 15 minutes.




Image of combi Impression coping and abutment assembly



Snap-on the plastic impression coping with the same sized diameter abutment



Impression taking  
Injection of impression material



Impression taking



Inner-surface of impression



Fabrication of provisional restoration or  
insertion of comfort cap

## LabSide



Seating of Lab analog



Confirm analog is secured in snap cap



Soft tissue model



Fabrication of master cast



Placement of burn-out cylinder



Consider the distance of opposing  
teeth, modify burn-out cylinders to its  
proper height.

# Abutment Level- Combi Abutment

[Multiple Units]



Connect the plastic bar in the middle of the trimmed burn-out cylinders to help support the resin pattern. Wax pattern may have shrinkage.



Wax-up



Completed framework



Trimming the extended margin with a rubber wheel



Metal framework and reamer



Removal of lip remnant with reamer caused by 'snap-on' mechanism



Metal Framework after removal of "Lip"



Metal coping adaptation (Completed framework)



Porcelain build-up final prosthesis

## Chairside




Insertion of final prosthesis and occlusal adjustment

\* If the combi analog is trimmed due to limited inter-occlusal space in the lab, make a reduction jig. Then a slight modify of the abutment in the oral cavity may be necessary to the height of the jig.

# Prosthetic Procedure 2

Impression Technique and Restoration Selection

## Dual / Milling / Angled / Metal-Casting / Temporary (Plastic & Ti) Abutment




# Fixture Level [Pick-up Type]- Dual Abutment

[Multiple Units]

## Clinical Procedure



## Chairside



# Fixture Level [Pick-up Type]- Dual Abutment

[Multiple Units]

## Laboratory Procedure



## Labside





# Fixture Level [Pick-up Type]- Dual Abutment

[Multiple Units]

## Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N-cm. Retighten after 15 minutes.



Insertion of the final prosthesis and occlusal adjustment

\* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.

## SCR- Labside



Formation of access hole with long transfer coping screw



Wax-up



Metal framework

## SCR- Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N-cm. Retighten after 15 minutes.



Insertion of final prosthesis and adjustment of occlusion

\* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.




# Fixture Level [Transfer Type]- Dual Abutment

[Multiple Units]

## Clinical Procedure



## Chairside



# Fixture Level [Transfer Type]- Dual Abutment

[Multiple Units]

## Laboratory Procedure



## Labside



Impression coping and analog connection. And insert impression coping into the impression.



Make sure the impression coping is fully seated into the impression



Soft tissue model



Fabrication of master cast



Soft tissue condition after the of impression coping



Measuring gingival height with depth gauge



Selection of Dual Abutment of proper diameter and gingival height



Verify by surveying the selected abutment. (Milling of the abutment is possible if necessary)



Fabrication of positioning jig

# Fixture Level [Transfer Type]- Dual Abutment

[Multiple Units]



Seat the cap with pattern resin



Completion of wax-up



Completion of metal framework

## Chairside



Final prosthesis built up on the framework with porcelain



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N-cm. Retighten after 15 minutes.



Insertion of final prosthesis, adjust occlusion place lab wax into opening of abutment to protect screw head then cement.

## SCR- Labside



Make an access hole in the resin cap by using the long open tray transfer screw.



Completed wax-up



Metal framework

## SCR- Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N-cm. Retighten after 15 minutes.



Insertion of final prosthesis and occlusal adjustment. Place wax into opening of the abutment prior to sealing with composite.

\* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply occlusal load on the prosthesis for 10~15 minutes.


# Fixture Level [Transfer Type]- Dual Milling Abutment

[Single Unit]

## Clinical Procedure



## Chairside



## Laboratory Procedure



# Fixture Level [Transfer Type]- Dual Milling Abutment

[Single Unit]

## Labside




Impression coping and analog connection. And insert impression coping into the impression.



Soft tissue model



Master cast



Selection of appropriate Dual Milling Abutment



Abutment after milling process



Fabrication of positioning jig



Fabrication of pattern resin cap



Completion of wax-up



Metal framework

## Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25~30N·cm. Retighten after 15 minutes.



Insertion of final prosthesis and occlusal adjustment


\* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In this case it is advised to apply acclusal load on the prosthesis for 10~15 minutes.




# Fixture Level [Pick-up Type]- Angled Abutment

[Single Unit]


## Clinical Procedure



## Chairside



## Laboratory Procedure





# Fixture Level [Pick-up Type]- Angled Abutment

[Single Unit]


## Labside



Impression coping with analog connections



Soft tissue formation and fabrication of master model



Unscrew then separate impression from the model.



Master cast



Select an Angled Abutment.



Modification of Angled Abutment & fabrication of positioning jig



Fabrication of pattern resin cap



Wax-up



Metal or zirconia framework

## Chairside



Final prosthesis



Insertion of the Angled Abutment using positioning jig



Insertion of final prosthesis and occlusal adjustment



***Super*** Líne


# Fixture Level- Metal-Casting Abutment

[Single Unit]

## Laboratory Procedure



## Labside




## Chairside



# Fixture Level [Pick-up Type]- Temporary Abutment

[Single Unit]



<Using Ti Abutment>



<Using Plastic Abutment>




Considering the opposing teeth before seating the Temporary Abutment, trim off the abutment as needed and complete the Temporary Abutment prosthesis with direct resin.

# Prosthetic Procedure 3

Impression Technique and Restoration Selection

## Screw Abutment

### Abutment Level Impression







# Abutment Level [Transfer Type]- Screw Abutment

[Multiple Units]

## Clinical Procedure



## Chairside






# Abutment Level [Transfer Type]- Screw Abutment

[Multiple Units]

## Laboratory Procedure



## Labside




# Abutment Level [Transfer Type]- Screw Abutment

[Multiple Units]



Metal framework



Removal of lip remnant in the interior of metal framework by using reamer



Completion of metal framework



Completion of porcelain



Insertion of final prosthesis and occlusal adjustment. Tighten with ti-retaining screw (10N-cm).

# Cementation Repair Method (SCRP)

[Screw & Cement Retained Prosthesis]

## In Light of Implant Prosthesis:

- A screw type restoration helps to simplify prosthesis repair, including insertion and removal of the prosthesis if necessary.
- Cement type restoration tend to have a stable occlusion and may enhance the adaptability. However the weak point is that it cannot be removed after permanent cementation.
- A Dual Abutment can be cemented or screw retained.
- Combi Abutments are cement retained and no occlusal hole is necessary.

## In Case of Screw Loosening or when Prosthesis Repair is Needed



In case of the following:  
screw loosening Prosthesis repair



In order to unscrew, form access hole on  
the occlusal surface using bur.



Unscrew, then remove the prosthesis  
from the oral cavity.



Both cemented prosthesis and  
abutment are removed.



Finish the repair then seat it inside the  
oral cavity.



Tighten the prosthesis with  
25~30N-cm by a screw driver.  
\* It is recommended that the abutment screw is  
retightened after 15 minutes.



Fill the access hole with cotton.



Fill the access hole with resin.



Final prosthesis

# Cementation Repair Method (SCRP)

[Screw & Cement Retained Prosthesis]

## Prosthesis Separation from Abutment due to Cement Loss



Remove the screw completely with screw driver and remove prosthesis from the patient's mouth.



Apply cement to the prosthesis.



Place it back into the patient's mouth.



After the cement setting, unscrew and remove the excessive cement.



Finish the repair and seat it inside the patient's mouth.



Tighten the prosthesis with 25~30N-cm with a screw driver.

## Adding to the Interproximal Contact Surface due to Prosthesis Loosening



Prosthesis loosening due to contact loosening



Form access hole using bur



Unscrew, then remove the cemented prosthesis with abutment in the oral cavity.




Contact adding with resin on the prepared under space.



Insert the prosthesis in the oral cavity and screw it in. Afterwards, perform light curing, then polish the contact area.



Position the prosthesis in the mouth and tighten the screw with 25~30N-cm, then fill up the access hole.




\* It is recommended that the abutment screw is retightened after 15 minutes.

# Prosthetic Procedure 4

Impression Technique and Restoration Type

## Overdenture Procedure

### Positioner / Mini Ball / Magnetic Attachment




# Positioner


## Chairside



Connect the Positioner Abutment onto the fixture.



Affix the impression coping on the Positioner Abutment.



Take Impression for the production of individual tray.




Produce the individual tray for denture impression.



After connecting the Positioner Abutment and the impression coping together, apply the impression material.



Take the final impression with the prepared individual tray.



After the impression material is set, discard the individual tray.





Image of the set final impression (with impression coping)

## Labside



Positioner Analog



Insert the Positioner Analog into the embedded impression coping.



Create the master model.



"Block out" procedure to achieve the space required for the metal socket.



Fabrication of denture with conventional method



# Positioner

## Case 1

## Chairside



Secure spaces for the female sockets.



Place the "block out spacer" on the Positioner Abutment in the patient's mouth.



Connect the metal socket onto the Positioner Abutment.



Apply a small amount of resin into the space created for the metal socket.



Position the denture in the mouth and wait until the resin is completely set.



Remove the white plastic socket (100gf) using the positioner tool and assemble with the regular plastic socket giving the desired retention force (300, 500 or 1000gf).



Remove the denture after the resin is fully set. Image of the denture with the metal socket.



Remove the block out spacer from the patient's mouth.



After polishing, the overdenture is completed.

## Case 2

## Chairside



Create holes for the placement of the metal sockets.



Place the "block out spacer" on the Positioner Abutment in the intraoral.



Connect the metal socket onto the Positioner Abutment.



Examine the interference between inner surface of the holes and the female sockets.



Apply the resin into the holes and wait until it is completely set.



Remove the white plastic socket (100gf) using the Positioner tool and assemble with the regular plastic socket giving the desired retention force (300, 500 or 1000gf).



Apply additional resin around the metal socket where there is a shortage of resin.



Apply resin around the metal socket.



After polishing, the overdenture is completed.

# Ball Attachment

## Case 1

## Chairside



Secure spaces for the female sockets.



Connect the female sockets to the Mini Ball Abutments in the intraoral.



Apply small amount of the resin into the secured area.



Position the denture in the mouth and wait until the resin is completely set.



Female sockets are placed in the denture.



After polishing, the overdenture is completed.

## Case 2

## Chairside



Create holes for the placement of the female sockets.



Connect the female sockets to the Mini Ball Abutments in the intraoral.



Examine the interference between inner surface of the holes and the female sockets.



Apply the resin into the holes and wait until it is completely set.



Place the female sockets.



Apply resin around the female sockets.



After polishing, the overdenture is completed.

# Magnetic Attachment

## Chairside



After Healing Abutment removal



Connect implant keeper with fixture and tighten it with 25~30N-cm.



Implant keepers connected with the fixtures



Position the magnetic assay on the implant keeper.



Secure spaces for the magnetic assays.



Examine the interference between inner divot of the denture and the magnets.

## Case 1



Apply resin on the divot of the denture's inner surface.



Position the denture into the mouth and wait until the resin is completely set.



Magnetic assays are placed in the denture.



Apply some of resin around the magnetic assays.



After the resin is completely set, remove excess. After polishing, the overdenture is completed.

# Magnetic Attachment

## Case 2



Create holes for the placement of the magnets.



Examine the interference between inner surface of the holes and the magnets.



Position the denture in the mouth and apply small amount of resin into the hole.



Wait until the resin is completely set.



After setting, remove denture from the mouth.



Add the resin around the magnets.



After polishing, the overdenture is completed.

# DENTIUM LONG-TERM CLINICAL DATA

2002

2003

2004

2005


2006

2007

2008



2002. 05. 17  
Pre-op



2002. 09. 04  
Post-op



2003. 03. 15  
Final prosthesis



# Dentium

2009

2010

2011

2012

2013

2014

2015

11 YEARS



2008. 04. 14  
5 years



2013. 12. 05  
11 years


Long term  
clinical cases data

Since

**2000**

OVER A **DECADE** OF  
COMMITMENT TO  
THE **BEST PRODUCTS**  
FOR DENTISTS AND  
PATIENTS





# New SuperLine II

## Product Catalog

**Dentium**

Specifications are subject to change without any notice.  
Some products listed in this catalog are not available in the market due to pending approval.

---

**HEAD OFFICE INDIA**

The Palm Springs, 8th Floor, Unit No 04, 05, 06, Palm Spring Plaza, Golf Course Rd, Sector 54, Gurugram, Haryana 122002

**HOMEPAGE**

[dentium.co.in](http://dentium.co.in)

**HEAD OFFICE**

174-10, #821 (Gangnam Ace Tower), Jagok-ro, Gangnam-gu, Seoul, Republic of Korea Tel +82-2-555-3750 Fax +82-2-6211-4681

**HOMEPAGE**

[www.dentium.co.kr](http://www.dentium.co.kr)