



bright Alone

Self-examination


dentium.co.in

Dentium

bright Alone


characteristic

- **Self-examination** can be conducted by using the **suction holder**.
- **Optimized moving distance** in limited space and **minimal A/S** with changed Dr. Table design
- **Convenient hygiene management** with detachable spittoon and instrument holder



Suction Holder

Suction holder enables easy **self-examination** and **auxiliary lighting** eliminates blind spots inside the mouth.




Detachable spittoon

Detachable spittoon enables easy hygiene management.



Floor Mount Type

Floor-mounted arm has **better durability** and enables the **efficient use of space**.



Option
Line-up

Table Type	Square Type
Table	400 x 275
Instrument	5 Holder High handpiece 2ea Low handpiece 1ea Scaler 1ea 3way syringe1ea





Option (C500)

Dental light with stable illumination and quality based
on **3-step color temperature control**
and **wide focusing pattern**

Foot Control_Function Type



Instrument

Line-up

Air Low Speed Handpiece

	FX22	FX65
Gear Ratio	1:1 Direct Drive	
For	CA Burs (Ø2.35mm)	HP Burs (Ø2.35mm)
Max Speed	30,000 min ⁻¹	40,000 min ⁻¹
Features	Clean Head System Push Button Chuck Ball Bearing	



High Speed Handpiece

M900L


Power	26 W
Speed	325,000~430,000 min ⁻¹
Spray	Quattro Spray
Features	Clean Head System Cellular Glass Optics Push Button Chuck

Scaler

Varios 170

Frequency	28~32 kHz
Max Output	11W
Dimensions	(W)63 x (D)54 x (H)26 mm
Module Weight	59 g





Chair Spec

Lift System	Hydraulic
Lift Height	490mm~890mm (400mm)
Recline range	0 ° ~70 ° (70 °)
Headrest	Double articulating
Dental light	LED light, Blue cut
Controller	Membrane Key + Color LCD
Spittoon	Manual rotation

BRAL_2021.04 [Rev.0]